

---

# **OL\_mrcnn**

***Release 1.0.0***

**Tom Soulaire**

**Jun 03, 2024**



# CONTENTS

<b>1</b>	<b>Installation</b>	<b>3</b>
<b>2</b>	<b>User Guide</b>	<b>5</b>
2.1	Image cropping . . . . .	5
2.2	Full pipeline . . . . .	5
2.3	Retraining your own model . . . . .	6
<b>3</b>	<b>Indices and tables</b>	<b>9</b>



This tool is a model based on a MRCNN architecture that enables to

- detect and crop cells in an image (grayscale, RGB, or more channels)
- classify cells in an image
- segment cells in an image



## INSTALLATION

To install the model on a local machine (requiring Python 3.9):

1. create conda environment:

```
conda env create -f environment.yml
```

2. activate conda environment:

```
conda activate OL_mrcnn
```

3. download the model weights [here](#) and place them in the folder /logs. These contain the original weights (trained on COCO dataset) and the weights trained on a custom dataset.





## 2.1 Image cropping

- Add your dataset in the folder `/data`
- OPTIONAL: preprocess your data with the `preprocessing.ipynb` notebook
- Configure the `image_cropper.ipynb` notebook:
  - `DEVICE`: device to use for inference. Default value is `'cpu:0'`.
  - `detection_min_confidence`: minimum confidence level for the detections. Default value is 0.7.
  - `detection_nms_threshold`: non-maximum suppression threshold. Eliminates the least confident detection when the IoU of 2 masks is above this value. Default value is 0.3.
  - `weights_subpath`: subpath in the `/logs` folder to the weights file.
  - `results_name`: name of the folder where the results will be saved.
  - `test_dir`: name of the folder where the images are stored.
  - `num_gpu`: number of GPUs to use for inference. Default value is 1.
  - `num_img_per_gpu`: number of images to process in parallel on each GPU. Default value is 1.
  - `VISUALIZE`: if True, displays the images with the detections. Default value is False.
- Run the notebook. The results will be saved in the folder `/results/results_name`.

## 2.2 Full pipeline

- In this setup, we run a first model to crop and classify objects in the images. Then we run a second model on the cropped images to get a refined mask.
- In the `model_pipeline.ipynb` notebook, configure the following parameters:
  - `DEVICE`: device to use for inference. Default value is `'cpu:0'`.
  - `gpu_count_macro`: number of GPUs to use for the first model. Default value is 1.
  - `num_img_per_gpu_macro`: number of images to process in parallel on each GPU for the first model. Default value is 1.
  - `min_confidence_macro`: minimum confidence level for the detections in the first model. Default value is 0.7.
  - `nms_threshold_macro`: non-maximum suppression threshold for the first model. Default value is 0.3.

- `nms_multiclass_macro`: non-maximum suppression threshold between classes for the first model. Default value is 0.3.
  - `gpu_count_micro`: number of GPUs to use for the second model. Default value is 1.
  - `num_img_per_gpu_micro`: number of images to process in parallel on each GPU for the second model. Default value is 1.
  - `min_confidence_micro`: minimum confidence level for the detections in the second model. Default value is 0.7.
  - `nms_threshold_micro`: non-maximum suppression threshold for the second model. Default value is 0.3.
  - `MACRO_MODEL_SUBPATH`: subpath in the `/logs` folder to the weights file of the first model.
  - `MICRO_MODEL_SUBPATH`: subpath in the `/logs` folder to the weights file of the second model.
  - `RESULTS_NAME`: name of the folder where the results will be saved.
  - `TEST_DIR`: name of the folder where the images are stored.
  - `VISUALIZE`: if True, displays the images with the detections. Default value is False.
- Run the notebook. The results will be saved in the folder `/results/RESULTS_NAME`.

## 2.3 Retraining your own model

### 2.3.1 Data structure

- Create a `/data` folder in the root directory.
- Inside the `/data` directory, put your images in a folder named `/imgs` and your binary masks in a folder named `/masks`. The name, size and format of the masks must match the images.
- In the `roi_labels_to_json.py` script, configure the `dir_path` in the `main()` function. Run in a terminal:

```
python roi_labels_to_json.py
```

- Move the label files to a `jsons` folder in the `data` directory.
- In the `format_data.py` script, configure the `dir_path` in the `main()` function. Configure the size of the training / validation / test datasets (usually 0.6, 0.2, 0.2). Run in a terminal:

```
python format_data.py
```

### 2.3.2 Retraining a single class model

- In the `custom.py` script, configure the following:
  - `GRAYSCALE`: if True, the model will be trained on grayscale images. Default value is False.
  - `DATA_PATH`: path to the dataset. Default value is `'/data'`.
  - `NAME`: name of the model.
  - `GPU_COUNT`: number of GPUs to use. Default value is 1.
  - `IMAGES_PER_GPU`: number of images to process in parallel on each GPU. Default value is 1.
  - `NUM_CLASSES`: number of classes. Default value is 2.

- EPOCHS: number of epochs. Default value is 50.
- STEPS PER EPOCH: number of steps per epoch. Default value is 50.
- LEARNING\_RATE: learning rate. Default value is 0.001.
- LAYERS: layers to train. Default value is 'heads'.
- DETECTION\_MIN\_CONFIDENCE: minimum confidence level for the detections. Default value is 0.7.
- DEVICE: device to use for training. Default value is 'cpu:0'.
- MAX\_GT\_INSTANCES: maximum number of instances in the ground truth. Default value is 100.
- DETECTION\_MAX\_INSTANCES: maximum number of instances in the detections. Default value is 35.
- in the ``CustomDataset`` class, modify or add lines :

- Run the script in a terminal:

```
python custom.py
```

### 2.3.3 Retraining a multi-class model

- Same instructions as before but on the `custom_multi.py` script.
- Run the script in a terminal:

```
python custom_multi.py
```



## INDICES AND TABLES

- `genindex`
- `modindex`
- `search`